

# Femosecond LIBS experiments from polymer targets


Rácz Péter

NAPLIFE collaboration  
Wigner Research Centre for Physics



# Laser Induced Breakdown Spectroscopy (LIBS)

LIBS: is a chemical analysis technology that uses a short (fs- ns) laser pulse to create a micro-plasma on the sample and based on atomic emission spectroscopy to measure elemental composition.





- Sample preparation-free measurements
- It can be used on any material, whether solid, liquid, or gas, foil
- Broad elemental coverage, including lighter elements, such as H, He Li, C, N, O, Na, and Mg, and heavy metallic elements too

<https://appliedspectra.com/technology/libs.html>

# Applications of LIBS

## Field of applications:

- Food science application: Measurement essential (Mg, Ca, and K) and toxic elements (Pb, Hg, St, Mn)
- Study of geological samples: identification of minerals
- Investigation of biological samples
- Pollution Monitoring
- Industrial application: in chemical industry material identification during manufacturing processes, in semiconductor industry used in semiconductor wafer and coating characterization and for quality control.
- The most appropriate applications of LIBS are in the nuclear and chemical industry, where quantitative or qualitative remote analysis, without any physical contact with the sample, is preferred.
- Further perspective in space exploration (for example on Mars rover)



# Applied laser systems for LIBS

## ns pulsed laser:

### -Nd:YAG Laser System:

- wavelength: 1064 nm, 532 nm
- pulse duration: 5–10 ns
- typical pulse energy: 10-300 mJ

### -The ruby laser:

- wavelength: 693 nm
- pulse duration: 20 ns;




## fs pulsed laser:

### - Ti:Sapphire chirped-pulse amplifier (CPA):

- wavelength: 800 nm
- pulse duration: 30–100 fs
- typical pulse energy: 10-50 mJ



## Advantages of fs-LIBS:

- low ablation threshold,
- improved spatial resolution for 3D mapping applications,
- small ablated mass, and reduced sample damage




# Laser Ablation Process

- Approximate time scales of nanosecond and femtosecond energy absorption and laser ablation along with various processes happening during and after the laser pulse is given



- Continuum light emission
- < 0.1 ~ 0.3 μsec







- Emission of discrete atomic lines at later times (~ 1 μsec).

# Timescales


- $\tau < 0.1 \sim 0.3 \mu\text{sec}$
- Plasma temperature 10000-20000 K
- Continuum light emission (Recombinations, ion lines)

- $\tau \sim 1 \mu\text{sec}$
- Plasma temperature 1000-2000 K
- Emission of discrete atomic lines





# Instrumentation for LIBS




# Applied laser system

Femtosecond Ti:Sapphire chirped-pulse amplifier



## Coherent Hydra-25:

- Pulse energy: max. 25-30 mJ
- Pulse length min. 40 fs
- Central wavelength 795 nm
- Repetition rate: 10 Hz
- Max Peak Power ca. 1 TW
- Focused max. peak intensity approx.  $10^{18}$  W/cm<sup>2</sup>





# Experimental Setup




## Spectrometer:

- LTB Demon spectrometer (Double Echelle Monochromator)
- Detector : ICCD
- Wavelength range: 190-900 nm
- Spectral resolution: 2.5-12 pm
- Simultaneous inspection range: 3 nm

## Environment conditions:




In vacuum ( $10^{-5}$  mbar)

Ar gas  $\sim$ 2-5 mbar for higher level signal



# Signal optimisation experiment

LIBS measurements from copper target





# Polymer targets

1. type


**UDMA:TEGDMA mixture (3:1)**

**UDMA** (urethane dimethacrylate),  $C_{23}H_{38}N_2O_8$

Used in dental medicine




**TEGDMA** (triethylene glycol dimethacrylate),  $C_{14}H_{22}O_6$



2. type

**Nanocomposite: UDMA:TEGDMA mixture (3:1) + AU nanorods**

Size of AU nanorods: 85 nm x 25 nm, Plasmonic resonance to 795 nm



3. type: Deuterated samples





**UDMA:MMA-D mixture (3:1) MMA-D (Methyl methacrylate):  $C_5D_8O_2$**

$$\frac{\text{No. of D atoms}}{\text{No. of H atoms}} = 0,32$$

See details in Attila Bonyár presentation

# Measurements from deuterated sample


- Balmer  $\alpha$  lines: hydrogen 656.28 nm and deuterium 656.11 nm.
- UDMA:MMA-D mixture (3:1)




# Measurements from bare polymers and nanocomposite samples

- Preliminary results
- Nanocomposite: UDMA:TEGDMA mixture (3:1) + AU nanorods

Single measurement  
25 mj, Ar 2 mbar  $\tau_d=0,9 \mu\text{s}$





Single measurement  
25 mj, Ar 2 mbar  $\tau_d=0,8 \mu\text{s}$





- Preliminary results
- Bare Polymer: UDMA:TEGDMA mixture (3:1)
- Nanocomposite: UDMA:TEGDMA mixture (3:1) + AU nanorods



# Summary

- LIBS experiment from deuterated polymer samples
- We can detect deuterium with LIBS method
- Measurements from nanocomposite targets (UDMA:TEGDMA mixture (3:1) + AU nanorods)
- Deuterium signal also present in this type of targets
- Possible explanation: deuterium production
- Next step: increasing the signal

Thank you for your attention!

Norbert Kroó, Tamás Biró, László Csernai, Miklós Veres, Márk Aladi,  
Miklós Kedves, Archana Kumári, Judit Kámán, István Rigó

Wigner Research Centre for Physics



Attila Bonyár, Alexandra Borók

Department of Electronics Technology, BME

Gábor Galbács

Department of Inorganic and Analytical Chemistry,  
University of Szeged







Concentration

|                              |             |
|------------------------------|-------------|
| $1.9 \times 10^{12}$ pice/ml | 0.1236 m/m% |
|------------------------------|-------------|